

Large Scale Microwave Heated SPPS of the ACP (65–74) Fragment Using Biotage® Initiator+ Alstra™

This application note demonstrates a variable scale synthesis on the Biotage® Initiator+ Alstra™ using the ACP (65–74) fragment as a model peptide.

Introduction

The ability to synthesize peptides in a range of different scales is a requirement in many research laboratories. The Biotage® Initiator+ Alstra™ fully automated microwave peptide synthesizer has a scale range of 5 µmol up to 2 mmol using three different reactor vial sizes (5 mL, 10 mL and 30 mL).

The acyl carrier protein fragment, ACP (65–74), H-Val-Gln-Ala-Ala-Ile-Asp-Tyr-Ile-Asn-Gly-NH₂ (VQAAIDYING-NH₂) (**1**), is a well-known so-called difficult sequence and is commonly used to evaluate the performance of new synthesis reagents and instrumentation.

Here we demonstrate a variable scale synthesis on the Initiator+ Alstra using the ACP (65–74) fragment as a model peptide.

Experimental

Materials

All materials were obtained from commercial suppliers; Sigma-Aldrich (Diisopropyl carbodiimide (DIC), piperidine, trifluoroacetic acid (TFA), triisopropylsilane (TIS), formic acid and Rink amide AM polystyrene resin, 1.1 mmol/g), Iris Biotech GmbH (Fmoc-amino acids), Fisher Scientific (NMP, DMF, diethyl ether and acetonitrile) and Novabiochem, Merck Millipore (Oxyma Pure and Rink amide AM polystyrene resin, 0.79 mmol/g). Milli-Q (Merck Millipore) water was used for LC-MS analysis.

N^ε-9-fluorenylmethoxycarbonyl (Fmoc) amino acids contained the following side-chain protecting groups: Asn(Trt), Asp(OtBu), Gln(Trt), Tyr(tBu).

Peptide Synthesis and Analysis

N^ε-Fmoc deprotection was performed at room temperature (RT) in two stages by treating the resin with 20% piperidine/NMP for 3 min followed by 20% piperidine/NMP for 10 min. The resin was then washed with NMP (4 ×). The peptides were synthesized using *N*^ε-Fmoc amino acids (3.0 eq., 0.5 M for 0.25 mmol scale and 3.0 eq., 0.7 M for 2 mmol scale), employing one of the following coupling protocols:

- » 0.25 mmol scale - DIC (3.0 eq., 0.5 M), Oxyma Pure (3.0 eq., 0.5 M) in NMP
- » 2.0 mmol scale - DIC (3.0 eq., neat), Oxyma Pure (3.0 eq., 4.5 M) in NMP

All couplings used microwave irradiation and were performed at 75 °C. After each coupling step, the resin was washed with NMP (4 ×).

For the 2 mmol scale synthesis a pause step was programmed in the synthesis after the coupling of the penultimate Gln residue to enable a small sample of the resin to be removed and checked by ninhydrin test, as it is known to be problematic. The ninhydrin test was slightly positive (some blue beads were observed), and therefore an additional Gln coupling step was programmed at this point “on-the-fly” and the synthesis continued until completion.

After the synthesis of the peptide sequence and final *N*^ε-Fmoc deprotection was completed, the resin was successively washed with NMP (3 ×), DCM (2 ×) and dried thoroughly. The peptides were cleaved from the solid support by treatment with TFA-TIS-H₂O (95:2.5:2.5) for 2 hours. The resin was separated by filtration and the cleavage cocktail was collected. The peptide was isolated by precipitation with cold diethyl ether (3 ×). The peptide was dissolved in 20%






Figure 1. RP-HPLC chromatogram and ESI-MS of 2 mmol scale synthesis of ACP (65-74) (**1**).

acetonitrile/water (the acetonitrile was added first), washed with DCM (2 x) and freeze dried to afford the product as a white powder. Analytical HPLC was performed on an Agilent 1100. The peptide was analyzed on a Resolux 200 Å C18 column (4.5 µm, 150 × 2.1 mm) with a flow rate of 1.0 mL/min. The following solvent system was used: solvent A, water containing 0.1% formic acid; solvent B, acetonitrile containing 0.1% formic acid. The column was eluted using a linear gradient from 0% buffer B to 90% buffer B over 20 min. Identification was carried out by ESI-MS (Agilent Technologies 6120 Qudropole LC/MS).

Results & Discussion

The ACP (65–74) peptide (**1**) was successfully synthesized using microwave irradiation at different scales with a crude purity of 97% (0.25 mmol scale) and 91% (2.0 mmol scale) (Figure 1) and confirmed by ESI-MS (Figure 1), calculated average isotopic composition for C₄₇H₇₅N₁₃O₁₅, 1061.55 Da. Found: m/z 1062.50 [M+H].

Conclusion

The ACP (65–74) peptide (**1**) was successfully synthesized in 2 mmol scale on a Biotage® Initiator+ Alstra™ fully automated microwave peptide synthesizer using standard pre-installed methods to afford the desired peptide in excellent crude purity.

We have demonstrated the advantages of using a robot liquid handler to dispense concentrated reagents which allows a lower excess of activated amino acids to be used, which is problematic for systems using nitrogen assisted valve based dispensing of reagents. With increasingly viscous reaction mixtures, mixing is an even more important factor at elevated temperatures for homogeneous heat distribution especially at larger scales, but this was overcome by the use of the highly effective oscillating mixer which enabled 5 min coupling times at 75 °C. In addition, the very flexible and intuitive software platform allowed various synthesis scales to be programmed effortlessly and the changing of a method “on-the-fly”.

EUROPE

Main Office: +46 18 565900
Toll Free: +800 18 565710
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com

NORTH AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
Order Fax: +1 434 296 8217
ordermailbox@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com

CHINA

Tel: +86 21 2898 6655
Fax: +86 21 2898 6153
cn_order@biotage.com

To locate a distributor,
please visit our website at
www.biotage.com